



**Caution: This Product Contains Natural Rubber Latex Which May Cause Allergic Reactions.**

The ShowerSafe cover is intended for use when bathing or showering. Do not submerge completely, or use when swimming or for recreational use. Use caution during use to prevent slipping.


**Cover desired area completely.**



**Wrap strap around limb making sure strap equally covers the cuff and the skin.**



**Stretch strap around limb and secure to white button.**



**Shower or bathe.**



1. Remove ShowerSafe from package, stretch cuff and completely cover area to be kept dry.
2. Position cover on limb so that white plastic button is easy to reach.
3. Stretch WrapAround strap around limb until snug and secure to white button.

**Hint:** It may be helpful to have someone assist in securing the strap.

**Note:** To assure a watertight seal, the WrapAround strap must equally overlap both the cuff and the skin.

**Important:** Overtightening the WrapAround strap may cause a lack of circulation. Remove immediately if numbness or tingling occurs. Do not wear for extended periods.

4. Shower or bathe. Do not submerge cover. Showering in a seated position is recommended when wearing a foot or leg cover.
5. After use, dry the outside of the cover before removing to prevent accidental wetting of the protected area.
6. **To dry**, turn cover inside out and hang with open end down and prevent pooling of water inside cover.

### Tips and Suggestions

- + Under certain conditions condensation may form inside cover during use. Lightly wrap limb with a thin towel to absorb moisture, if desired.
- + If the ShowerSafe cover is punctured or torn, repair with waterproof tape such as duct tape or first aid tape.