

ShowerSafe is a reusable waterproof plastic covering that keeps casts, bandages and dressings dry during bathing and showering.

Patients with fractures, wounds and IV sites are frequently unable to bathe easily due to concerns about water damage to casts, bandages and dressings.

The ShowerSafe is made of heavy, 4 mil plastic for maximum **durability** and protection.

Unique **DoubleSeal** design protects in two ways- the first seal is made by the elastic cuff that adjusts to limb size and provides a snug fit, the second seal is made by the integral **WrapAround** strap that seals around the elastic cuff and the limb for the ultimate in waterproof protection.

The use of the ShowerSafe allows patients to *safely* enjoy a bath or shower with confidence:

- Prevents cast replacement
- Prevents bandage replacement
- Prevents dressing changes
- Prevent infection from moisture exposure
- Prevents risk of injury from trying to keep limbs outside the tub or shower

ShowerSafe is available in nine sizes, including our unique knee and elbow models to meet a complete range of applications including:

- Casts
- Bandages and Dressings
- Surgical Wounds
- I.V. Sites
- PICC Catheter Lines
- Orthopedic procedures



Order #	Size (width x length)	Fits
4400-S	Small Arm (8" x 30")	½ Adult arm, child up to 85 lbs
4400-M	Medium Arm (8" x 30")	Adult up to 150 lbs
4400-L	Large Arm (10" x 30")	Adult over 150 lbs
4410-S	Small Leg (8" x 30")	½ Adult leg, child up to 85 lbs
4410-M	Medium Leg (12" x 43")	Adult up to 170 lbs
4410-L	Large Leg (15" x 43")	Adult over 170 lbs
4430-S	Small Knee & Elbow (9" x 16"; 5" opening)	Adult up to 150 lbs
4430-L	Large Knee & Elbow (9" x 24"; 5" opening)	Adult over 150 lbs
4420	Foot and Ankle (10" x 20")	Universal size fits to just below knee

An anti-itch solution for relief of cast discomfort.

Castblast is a talc spray which instantly and safely relieves the itching and irritation associated with casting.



Order #	Description
4500	Castblast Talc Spray, 4.5 oz. bottle