

Product Description	Sterile Catheters Order #	Non-Sterile Catheters Order #
Res-Q-Vac Hand Powered Pump, Reusable	—	2214
Res-Q-Vac Combination Kit: Res-Q-Vac Pump, Canister w/ Rigid Yankauer Catheter <i>Replacement Canister/Catheter Kit (case of 10)</i>	2242 2215	2248 2251
Res-Q-Vac Combination Kit: Res-Q-Vac Pump, Canister, Soft Large Bore Catheter <i>Replacement Canister/Catheter Kit (case of 10)</i>	2243 2216	*
Res-Q-Vac Combination Kit: Res-Q-Vac Pump, Canister, Soft Large Bore Catheter, and 3.5" and 14" 10 Fr Catheters <i>Replacement Canister/Catheter Kit (case of 10)</i>	2221 2231	2240 2250
Res-Q-Vac Combination Kit: Res-Q-Vac Pump, Canister, 3.5" and 14" 14 Fr Catheters <i>Replacement Canister/Catheter Kit (case of 10)</i>	2222 2223	*
Res-Q-Vac Combination Kit: Res-Q-Vac Pump, Canister, 3.5" and 14" 10 Fr Catheters <i>Replacement Canister/Catheter Kit (case of 10)</i>	*	2253 2255
Product Specifications <ul style="list-style-type: none"> • Res-Q-Vac pump and canister are non-sterile • Canister Size: 300 cc • Soft Large bore catheter: 1/2" (38 Fr) diameter by 7" long • Rigid Yankauer catheter: 5/16" (24 Fr) diameter by 7" long 	* Other combinations may be available by special order. Call Customer Service at 800-325-9044	