

# ShowerSafe™

# Cast & Bandage Cover

## Sizing Chart

(sizes are approximate)

	Arm	Leg	Knee & Elbow	Foot & Ankle
Small	1/2 arm adult, Full arm child to 85 lbs. 5" diameter x 21" long	1/2 leg adult, Full leg child to 85 lbs. 6.5" diameter x 26" long	Up to 150 lbs. 5" diameter x 14" long	Universal 6.5" diameter by 17" long
Medium	Full arm to 150 lbs. 5" diameter x 26" long	Full leg to 170 lbs. 7.5" diameter x 39" long	■ ■ ■ ■ ■	
Large	Full arm, adult over 150 lbs. 6.5" diameter x 26" long	Full leg, adult over 170 lbs. 9.5" diameter x 39" long	Over 150 lbs. 6.5" diameter x 14" long	

Caution: This Product Contains Natural Rubber Latex Which May Cause Allergic Reactions.

 **Trademark  
Medical**

Trademark Medical  
St. Louis, MO 63011  
800-325-9044  
[www.trademarkmedical.com](http://www.trademarkmedical.com)