

Res-Q-Vac is a hand-powered portable suction system that creates hospital-level suction power without batteries or electricity.

Res-Q-Vac is a hand-powered portable suction system consisting of the reusable Res-Q-Vac suction pump handle, and a variety of disposable suction canister and suction tube configurations. Res-Q-Vac is always ready to use on a moment's notice because it does not rely on batteries or electricity to work- and it requires no maintenance!

- The Res-Q-Vac weighs less than ½ pound and uses no batteries or electricity
- Two compressions of the handle creates hospital-level suction power of over 400 mmHg. Suction levels of over 550 mmHg can be achieved if needed.

Res-Q-Vac is ideal for oral suctioning, and other situations where suction is needed:

- Crash carts and emergency response kits, the emergency room, isolation units, for general floor use and as a back up suction system throughout the hospital.
- Dental use as a back-up suction system in the event of power or suction failure
- Gentle, controlled suction during childbirth and for meconium aspiration in newborns, and for pediatric emergencies.
- Nursing homes and in extended care settings where wall suction may not be available

Res-Q-Vac is available in sterile or non-sterile, or with FULL STOP PROTECTION.

Res-Q-Vac is also available in combination with our Plak-Vac Oral Suction Brush to allow effective oral hygiene for people with swallowing difficulties. (See photo below)



Res-Q-Vac with soft large bore catheter



Res-Q-Vac with rigid yankauer tube



Res-Q-Vac with home suction kit



Plak-Vac/
Res-Q-Vac
combination



Full Stop Protection

Acts like a .22 micron filter to prevent exposure to blood and other potentially infectious materials to help maintain compliance with OSHA Bloodborne Pathogens regulations.



Description	FULL STOP PROTECTION NON-STERILE Order #	STANDARD PRODUCT NON-STERILE Order #	STANDARD PRODUCT STERILE Order #
Res-Q-Vac Reusable Pump	–	2214	–
Res-Q-Vac Combination Kits (pump, canister, catheter, cap & label)			
Rigid Yankauer Kit	2242	2248	2220
Soft Large Bore 36 Fr. Catheter Kit	2243	2240	2249
Rigid Yankauer Kit with ET Adapter Kit	2244	–	–
Rigid Yankauer Kit with ET Adapter and 10 Fr. 14" Catheter Kit	–	2219	2241
Soft Large Bore Catheter Kit with 10 Fr. Catheter Kit	–	–	2221
Neonatal/Pediatric Kit with 10 Fr. 3.5" Catheter, and 10 Fr. 14" Catheter,	–	–	2253
Home Suction Kit with 14 Fr. 3.5" Catheter, 14 Fr. 14" Catheter, and Carry Pouch	–	2222	–
Replacement Canister Kits (canister, catheter, cap & label -10 units/case)			
Soft Large Bore 36 Fr. Catheter	2245	2250	2216
Rigid Yankauer Catheter	2246	2251	2215
Rigid Yankauer Catheter w/ET Adapter	2247	2252	2217
Neonatal/Pediatric with 10 Fr. 3.5" Catheter, and 10 Fr. 14" Catheter,	–	–	2254
Home Suction Kit with 14 Fr. 3.5" Catheter and 14 Fr. 14" Catheter	–	2223	–