


# Wound Measuring Guide

Wound Management Aid

The disposable wound measuring guide provides a convenient and economical way to document and monitor wound healing.

The single use, disposable measuring guide enables convenient and effective measurement of dermatological wound sites such as skin lesions, pressure sores, stomas, wounds and incisions while eliminating cross contamination. Accurate documentation of wound length, width and diameter is now made easy and economical with dual scale ruler templates. Compliant 2 mil plastic conforms to curved surfaces for greater accuracy. Packed in a convenient dispenser box.



Order #	Description
B250	Disposable Measuring Guide, 3-3/4" x 5-1/2" - 250 count
B500	Disposable Measuring Guide, 3-3/4" x 5-1/2" - 500 count
B10-250	Large Disposable Measuring Guide, 5-1/4" x 6-1/2" - 250 count

# Heel Protector

Foam Support Aid

Heel Protector surrounds the injured heel and ankle with soft and firm convoluted foam. Secured by straps the heel protector offers total support without limiting mobility. Two per package.



Order #	Description
4700	Heel Protector - High Density Foam - pair