

Res-Q-Vac is a hand-powered portable suction system that creates hospital-level suction power without batteries or electricity.

Res-Q-Vac is a hand-powered portable suction system consisting of the reusable Res-Q-Vac suction pump handle, and a variety of disposable suction canister and suction tube configurations. Res-Q-Vac is always ready to use on a moment's notice because it does not rely on batteries or electricity to work- and it requires no maintenance!

- The Res-Q-Vac weighs less than ½ pound and uses no batteries or electricity
- Two compressions of the handle creates hospital-level suction power of over 400 mmHg. Suction levels of over 550 mmHg can be achieved if needed.

Res-Q-Vac is ideal for oral suctioning, and other situations where suction is needed:

- Crash carts and emergency response kits, the emergency room, isolation units, for general floor use and as a back up suction system throughout the hospital.
- Dental use as a back-up suction system in the event of power or suction failure
- Gentle, controlled suction during childbirth and for meconium aspiration in newborns, and for pediatric emergencies.
- Nursing homes and in extended care settings where wall suction may not be available

Res-Q-Vac is available with sterile or non-sterile catheters and yankauers.



Res-Q-Vac with soft large bore catheter



Res-Q-Vac with rigid yankauer catheter



Res-Q-Vac with 3.5" & 14" 14Fr catheters



Res-Q-Vac Pump Handle

Res-Q-Vac is also available in combination with our Plak-Vac Suction Toothbrush to allow effective oral hygiene for people with swallowing difficulties.



Plak-Vac/Res-Q-Vac combination

Product Description	Sterile Catheters Order #	Non-Sterile Catheters Order #
Res-Q-Vac Hand Powered Pump, Reusable	—	2214
Res-Q-Vac Combination Kit: Res-Q-Vac Pump, Canister w/ Rigid Yankauer Catheter <i>Replacement Canister/Catheter Kit (case of 10)</i>	2242 2215	2248 2251
Res-Q-Vac Combination Kit: Res-Q-Vac Pump, Canister, Soft Large Bore Catheter <i>Replacement Canister/Catheter Kit (case of 10)</i>	2243 2216	*
Res-Q-Vac Combination Kit: Res-Q-Vac Pump, Canister, Soft Large Bore Catheter, and 3.5" and 14" 10 Fr Catheters <i>Replacement Canister/Catheter Kit (case of 10)</i>	2221 2231	2240 2250
Res-Q-Vac Combination Kit: Res-Q-Vac Pump, Canister, 3.5" and 14" 14 Fr Catheters <i>Replacement Canister/Catheter Kit (case of 10)</i>	2222 2223	*
Res-Q-Vac Combination Kit: Res-Q-Vac Pump, Canister, 3.5" and 14" 10 Fr Catheters <i>Replacement Canister/Catheter Kit (case of 10)</i>	* 	2253 2255
Product Specifications <ul style="list-style-type: none"> • Res-Q-Vac pump and canister are non-sterile • Canister Size: 300 cc • Soft Large bore catheter: 1/2" (38 Fr) diameter by 7" long • Rigid Yankauer catheter: 5/16" (24 Fr) diameter by 7" long 	* Other combinations may be available by special order. Call Customer Service at 800-325-9044	