

Face Fit Goggles

Protective
Eyewear

Standard and Anti-Fog enhanced face fit goggles for comfortable protection

Face fitting goggles which are comfortable and soft, fashioned from clear, shatter resistant plastic providing distortion-free vision. The eyes remain completely covered and protected from splatter, particulate and debris. Adjustable strap maintains secure fit and allows for wear over normal eyeglasses.



Order #	Description
9552	Anti-Fog Clear Face Fit Goggles
9553	Standard Clear Face Fit Goggles

Blade-Safe

Surgical
Blade Remover

The Blade-Safe is designed to prevent sharps injury when removing surgical blades from scalpel handles.

The Blade-Safe allows safe removal of scalpel blades in three steps:

- Place scalpel blade in Blade-Safe device
- Snap Blade-Safe lid closed to contain blade
- Remove blade by holding Blade-Safe in one hand and pulling scalpel handle with the other.

After handle has been removed, the scalpel blade is contained within the Blade-Safe. Dispose of the Blade-Safe and blade in appropriate sharps container.

Use for blades up to and including #20, and for #3 and #4 reusable handles only.

Available in single sterile and non-sterile packages.



Order #	Description
730	Blade-Safe Surgical Blade Remover - non-sterile individually packaged
750	Blade-Safe Surgical Blade Remover - sterile individually packaged