

The Sani-Stor® instrument holder provides a simple, convenient and hygienic way to store instruments between uses. Made of rigid translucent polypropylene, the Sani-Stor holder protects the patient by guarding the instrument against accidental contact with potentially contaminated surfaces such as bed linens, furniture and equipment, and protects staff and family members by preventing accidental contact with potentially contaminated instruments.

The Sani-Stor holder allows air to circulate around instruments during storage and between uses, promoting air drying and helping to prevent moist conditions which allow for bacterial growth.

The Sani-Stor holder is ideal for storing any patient care instrument that needs to be available for use at the point of care, particularly oral suction instruments, such as the yankauer.

Proper handling and storage of the yankauer is more important than ever as oral suction instruments have been shown to become colonized with pathogenic organisms that could lead to hospital acquired infections.\*

**Sani-Stor Features**

**Multi-Port Locking Slot**



Secures and keeps instruments above bottom of holder.

**Adjustable Elastomeric Band**

Secures holder to rail and resists rotating out of position.

**Translucent Polypropylene**

Permits monitoring of holder for any build-up of debris.



Order #	Description
2650	Sani-Stor Instrument Holder 50/case

\*Sole ML et al. Bacterial Growth in Secretions and in Suctioning Equipment of Orally Intubated Patients: A Pilot Study. Am J Crit Care. 2002;11:141-49.

**LATEX FREE**

**Covered by U.S. Patent No. 8,245,857**