

Plak-Vac® is a reusable single patient oral hygiene instrument for patients with difficulty swallowing and expectorating, or who cannot perform routine oral care.

The Plak-Vac Suction Toothbrush provides the ultimate in care for patients requiring assistance with oral hygiene. A combination suction instrument and toothbrush, the Plak-Vac allows the removal of debris, plaque, bacteria and fluid from the mouth of patients who are unable to swallow or have difficulty swallowing and expectorating. The Plak-Vac has soft bristles that are gentle to gums, teeth and other tissues in the mouth. Suction level is controlled by a finger control port located in the Plak-Vac handle.

Use the Plak-Vac to:

- Make oral hygiene simple and effective.
- Remove mucus and food build-up, as well as plaque, bacteria and fluid.
- Help eliminate chronic halitosis.
- Reduce the risk of infection, aspiration and pneumonia.
- Eliminate the need for separate yankauers, toothbrushes and suction catheters.

Plak-Vac is ideally suited for use with both wall suction and portable pumps.

Our Plak-Vac/Res-Q-Vac combination features the Plak-Vac suction toothbrush with the Res-Q-Vac hand powered suction system.

It is the perfect system for home care and long term care use where wall or portable suction may not be available, and since it does not require a power source it is completely portable.



Order #	Description
2200	Plak-Vac Suction Toothbrush
2213	Plak-Vac/Res-Q-Vac Combination
2201	Plak-Vac Suction Toothbrush for use with Res-Q-Vac Hand Operated Pump



The **Plak-Vac Storage Holder** provides a simple, convenient and hygienic way to store the Plak-Vac suction toothbrush between uses. The transparent Plak-Vac storage holder helps avoid contamination of the Plak-Vac toothbrush by preventing accidental contact with bed linens, furniture and other potentially contaminated surfaces- and permits air drying. The Plak-Vac storage holder is great for storing other patient care devices such as yankauer suction instruments, as well. Use with the Plak-Vac storage holder bracket or clip.



The **Plak-Vac Storage Holder Bracket** is a permanent metal receptacle for the Plak-Vac storage holder. Coated with a durable anti-microbial paint, the Plak-Vac storage holder bracket mounts permanently to walls, carts or any desired surface. Mounting hardware not included.



The **Plak-Vac Storage Holder Clip** offers a semi-permanent way to secure the Plak-Vac storage holder. Made of durable injection molded plastic, the Plak-Vac storage holder clip mounts to smooth surfaces with removable two-sided tape for use where a permanent installation is not desirable or feasible.



Plak-Vac Oral Care Kit

Trademark Medical also offers a customized Plak-Vac Oral Care Kit Program providing a dedicated component formulary to assure the implementation of your protocol while minimizing cost and waste. Contact Trademark Medical for more details.

Order #	Description
2700	Plak-Vac Storage Holder
2701	Plak-Vac Storage Holder Bracket
2702	Plak-Vac Storage Holder Clip